

# SSIbio

## Bullet Boxes

### New! 5 mL Tube Storage & Freezer Box With Lid

SSI is pleased to introduce a new 5 mL tube rack, expanding its line of handy storage solutions for most tubes in the lab.

- Three models accommodate common tube sizes (0.5, 1.5/2.0, and 5.0 mL).
- Ideal for screw cap microtubes and microcentrifuge tubes.
- Small footprint saves bench and freezer space.
- Side notches facilitate easy lid removal.
- Suitable for temperatures as low as -90°C as well as autoclaving.
- Racks come with lid.

### Dimensions

- 5 mL Box: 143.8 x 75.2 x 53.0 mm
- 1.5/2.0 mL Box: 143.8 x 75.2 x 53.0 mm
- 0.5 mL Box: 139.7 x 71.1 x 61.8 mm

Item #	Tube Volume	Description	Color	Packaging
5527-00	5 mL	21-Place Freezer Box	Natural	5 racks/unit, 4 units/case
5520-00	1.5 – 2.0 mL	50-Place Freezer Box	Natural	
5530-00	0.5 mL	50-Place Freezer Box	Natural	

