

Vertex Macro Volume Tip Compatibility Chart

- Improved accuracy over serological pipettes when measuring large volumes of solution.
- Tip length facilitates access into bottles, cell culture flasks and tubes.
- Free of detectable RNase, DNase, DNA, and PCR inhibitors. Sterile versions are also pyrogen-free.



Macro Tips Description

Item #	Volume	Description	Packaging
4410-00	5 mL	Gilson-fit, Bulk, Natural	250/unit, 10 units/case
4410-0F	5 mL	Gilson-fit, Bulk, Filtered, Natural	100/unit, 10 units/case
4420-00	5 mL	Eppendorf-fit, Bulk, Natural	250/unit, 10 units/case
4420-0F	5 mL	Eppendorf-fit, Bulk, Filtered, Natural	100/unit, 10 units/case
4411-00	5 mL	Gilson-fit, Racked, Natural	50 tips/unit, 10 units/case
4421-00	5 mL	Eppendorf-fit, Racked, Natural	50 tips/unit, 10 units/case
4411-S0	5 mL	Gilson-fit, Racked, Sterile, Natural	50 tips/unit, 10 units/case
4421-S0	5 mL	Eppendorf-fit, Racked, Sterile, Natural	50 tips/unit, 10 units/case
4411-SF	5 mL	Gilson-fit, Racked, Sterile, Filtered, Natural	50 tips/unit, 10 units/case
4421-SF	5 mL	Eppendorf-fit, Racked, Sterile, Filtered, Natural	50 tips/unit, 10 units/case
<hr/>			
4520-01	10 mL	Bulk, Natural	200/unit, 5 units/case
4530-00	10 mL	Gilson-fit, Bulk, Natural	200/unit, 5 units/case
4530-0F	10 mL	Gilson-fit, Bulk, Filtered, Natural	50/unit, 10 units/case
4521-01	10 mL	Racked, Natural	25 tips/rack, 10 racks/case
4531-00	10 mL	Gilson-fit, Racked, Natural	25 tips/rack, 10 racks/case
4521-S1	10 mL	Racked, Sterile, Natural	25 tips/rack, 10 racks/case
4531-S0	10 mL	Gilson-fit, Racked, Sterile, Natural	25 tips/rack, 10 racks/case
4521-SF	10 mL	Racked, Sterile, Filtered, Natural	25 tips/rack, 10 racks/case
4531-SF	10 mL	Gilson-fit, Racked, Sterile, Filtered, Natural	25 tips/rack, 10 racks/case



Macro Tips Compatibility

Item #	Eppendorf	Socorex	Gilson	Rainin	Accumax	Biohit	Finnpipette
4410-00			✔	✔			
4410-0F			✔	✔			
4420-00	✔	✔				✔*	
4420-0F	✔	✔					
4411-00			✔	✔			
4421-00	✔	✔				✔*	
4411-S0			✔	✔			
4421-S0	✔	✔				✔*	
4411-SF			✔	✔			
4421-SF	✔	✔				✔*	
<hr/>							
4520-01		✔			✔		✔
4530-00	✔		✔	✔		✔*	
4530-0F	✔		✔	✔		✔*	
4521-01		✔			✔		✔
4531-00	✔		✔	✔		✔*	
4521-S1		✔			✔		✔
4531-S0	✔		✔	✔		✔*	
4521-SF		✔			✔		✔
4531-SF	✔		✔	✔		✔*	

*Does not fit Proline pipettors, but does fit all other Biohit pipettors, including Proline Plus.