




**NOTE(S): UNLESS OTHERWISE SPECIFIED.**

1. MATERIAL: 1.4mil ALUMINUM FOIL FILM WITH 1.1mil ACRYLIC ADHESIVE ON A 4.6mil BROWN POLY-COATED PAPER RELEASE LINER.
2. BACKSCORES THROUGH 4.6mil BROWN POLY-COATED PAPER RELEASE LINER ONLY.
3. PERFORATIONS THROUGH 1.4mil ALUMINUM FOIL FILM ONLY.



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS	CAD GENERATED DRAWING, DO NOT MANUALLY UPDATE		 <b>1310 Thurman Street</b> <b>Lodi CA, 95240</b> <b>Ph: (209) 333-2120</b> <b>Fx: (209) 333-8623</b> <b>info@scientificspecialties.com</b>
	APPROVALS	DATE	
SSI MATERIAL:	DRAWN: P. DO	2/22/2018	<b>3640-00-XL</b> <b>ALUMASEAL II</b>
MANUFACTURER'S MATERIAL:	CHECKED:		
FINISH:	DESIGN ENG:		SIZE: SHRINK RATE:
DO NOT SCALE DRAWING	MFG ENG:		<b>A</b>
	QUAL ENG:		REV.: REL
			SCALE: 1:2
			JOB NO.:
			SHEET: 1 OF 1

PRINTED 2/22/2018