


ITEM	PART NUMBER	DESCRIPTION	QTY.
1	6121-XX	1.8mL INT. THD. CRYOVIAL TUBE	1
2	M9-0-2285	SILICONE CAP SEAL	1
3	6100-XX	EXT. THD. CRYOVIAL CAP	1

CAD GENERATED DRAWING. DO NOT MANUALLY UPDATE		 1310 Thurman Street Lodi CA, 95240 Ph: (209) 333-2120 Fx: (209) 333-8623 info@scientificspecialties.com				
APPROVALS	DATE					
DRAWN S.Gallagher	9/1/11	<h1>6121-S0</h1> <h2>1.8mL CRYOVIAL ASSEMBLY</h2>				
CHECKED						
RESP ENG						
MFG ENG						
FINISH	---	QUAL ENG	QUAL ENG	SIZE A	SHRINK RATE	REV. REL
DO NOT SCALE DRAWING		SCALE	JOB NO.	SHEET	1 OF 1	