



SECTION A-A



ITEM	PART NUMBER	DESCRIPTION	QTY.
1	6122-XX	1.8mL INT. THD. CRYOVIAL TUBE	1
2	M9-0-2285	SILICONE CAP SEAL	1
3	6100-XX	EXT. THD. CRYOVIAL CAP	1

CAD GENERATED DRAWING. DO NOT MANUALLY UPDATE		 1310 Thurman Street Lodi CA, 95240 Ph: (209) 333-2120 Fx: (209) 333-8623 info@scientificspecialties.com			
APPROVALS	DATE				
DRAWN S.Gallagher	9/1/11	<h1 style="text-align: center;">6122-S0</h1> <h2 style="text-align: center;">1.8mL CRYOVIAL ASSEMBLY</h2>			
CHECKED					
SSi MATERIAL: ---	RESP ENG				
MANUFACTURERS MATERIAL: ---	MFG ENG				
FINISH ---	QUAL ENG	<table border="1"> <tr> <td>SIZE A</td> <td>SHRINK RATE</td> <td>REV. REL</td> </tr> </table>	SIZE A	SHRINK RATE	REV. REL
SIZE A	SHRINK RATE	REV. REL			
DO NOT SCALE DRAWING		<table border="1"> <tr> <td>SCALE</td> <td>JOB NO.</td> <td>SHEET 1 OF 1</td> </tr> </table>	SCALE	JOB NO.	SHEET 1 OF 1
SCALE	JOB NO.	SHEET 1 OF 1			

PRINTED 9/2/2011