







SECTION A-A



ITEM	PART NUMBER	DESCRIPTION	QTY.
1	6112-XX	1.0mL INT. THD. CRYOVIAL TUBE	1
2	M9-0-2285	SILICONE CAP SEAL	1
3	6100-XX	EXT. THD. CRYOVIAL CAP	1

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:

FRACTIONS	DECIMALS	ANGLES
$\pm 1/64$	$.XX \pm .01$	± 0.1
	$.XXX \pm .001$	
	$.XXXX \pm .0005$	

SSI MATERIAL: ---

MANUFACTURERS MATERIAL: ---

FINISH ---

DO NOT SCALE DRAWING

CAD GENERATED DRAWING,
DO NOT MANUALLY UPDATE

APPROVALS	DATE
DRAWN S. Gallagher	9/1/11
CHECKED	
RESP ENG	
MFG ENG	
QUAL ENG	



1310 Thurman Street
Lodi CA, 95240
Ph: (209) 333-2120
F: (209) 333-8623
info@scientificspecialties.com

6112-S0
1.0mL CRYOVIAL ASSEMBLY

SIZE	SHRINK RATE	REV.
A		REL
SCALE	JOB NO.	SHEET 1 OF 1